


CrystalWave™

LUMETRICS® offers CrystalWave™ dry testing of intra-ocular lenses. Focal length, wavefront, and MTF measurements are combined in a single instrument using a Shack-Hartmann wavefront sensor.

For normal software operations, a simple user interface guides the alignment of the DL and the measurements. Researchers can access the CLAS-2D™ software for advanced data analysis, including Zernike polynomials, coma, and spherical aberration.

The CrystalWave also supports production line installations and is insensitive to vibrations. The measurement zone accepts a single DL or multiple DLs in a carousel. The optional CrystalServer™ software allows user written programs to control CrystalWave in coordination with robotic equipment.

The combined measurement of power, wavefront, and MTF measurements in one instrument minimizes part handling. Older technologies required separate instruments for power and MTF and required subjective judgments by technicians. Now wavefront technology provides objective measurements that are repeatable for process control.



A unique algorithm removes the usual dependence of measured power on the axial location of the DL. This algorithm eliminates the need for mechanical restraints on the DL and eliminates the resulting distortion of soft materials. CrystalWave also includes an inspection camera with back illumination of the DL for monitoring scratches, tears, irregular edges, debris, and haptics orientation.





LUMETRICS®

Precision Measurement Solutions

Specifications:

CrystalWave Sensor:

101x101 lenslet array grid
 Field of View : 7.4 x 7.4 mm
 Spatial Resolution : 0.072 mm
 Spherical Range : 7 to 120 Diopters (dry only version)
 Cylinder Range : 6 diopters for 6 mm analysis diameter
 XY translation stages for adjustment of DL (no cuvette)

Standard Equipment:

Inspection Camera displayed on controller monitor
 DL illumination LED ring
 Rotating Dry Mount for DL having monofilament haptics
 Controller with monitor
 Configured for vertical operation

Mechanical:

Height, width, length : 35 x 11 x 13 inches
 CrystalWave™ may be used vertically only

Accuracy and Repeatability:

Power Accuracy : 0.15 Diopters Sphere and Cylinder (Range 6-40D)
 Power Repeatability : 0.05 Diopters
 Sphere and Cylinder Spherical Aberration Accuracy : 0.004 microns
 (Z42) Spherical Aberration Repeatability : 0.002 microns (Z42)

Data Outputs:

Sphere/Cylinder/Axis
 Focal Length
 Spherical Aberration
 MTF of DL itself (calculated from wavefront)
 MTF of DL installed in ISO mode eye (optional) (calculated from wavefront)
 Zernike Polynomials
 Zonal Wavefront Map
 Point Spread Function
 Power Map

The screenshot shows the Crystalwave software interface with the following data:

- Lot #: QA7H9
- Lens ID #: 12
- Operator: JR
- Date / Time: T 10/20/2004 15:20:27.187000
- Media: Air, Water
- Collimation Error: -0.037 Dpt
- Astigmatism: 0.056 Dpt
- XAng: -10.148 mr
- YAng: -4.394 mr
- Z22: 0.03, Z31: 0.04, Z32: 0.01
- Analysis Diameter: 5.0 mm
- PV WFE: 0.3 microns
- Spherical Aberration: -0.014 microns
- Power: 11.173 D
- File Name: QA7H9_12_102020041520270011
- Cuvette Temperature: NA C
- Mount Temperature: NA C

Sample Screen Shot

Options:

CrystalServerActive-X automation DLL
 callable from Labview, VisualBasic, etc

Specifications are subject to change without notice.
 CLAS-2D, ClearWave, and CrystalWave are trademarks of FOM WaveFront Sciences, LLC